

Fisher Lab Apoptosis Fluorescent Protocol

Solutions

1. Make working strength TDT enzyme. Mix **77ul** of reaction buffer with **33ul** of TDT enzyme in **390ul** of PBS. Mix well by vortexing.
2. Make Stop-Wash buffer. (**1ml** of Stop Wash buffer into **34ml** of distilled water.)

Procedure (Keep wrapped in foil throughout the procedure)

1. Thoroughly wash specimen in PBS (minimum of 2 hrs. with several washes in between.)
2. Immediately apply **75ul** equilibration buffer in addition to **425 ul** of PBS for a total of **500 ul** of solution directly onto the specimen. Incubate for 3 minutes at room temperature.
3. Remove equilibration buffer, apply working strength TDT enzyme (solution #1), wrap in aluminum foil and label (Fisher Lab x 3611, Name, Date) and place in incubator at 37 degree (located on 5th floor) for 1 hour.
4. Remove TDT enzyme, and apply stop-wash buffer (solution #2). Gently shake for 30 seconds, followed by 10 min. incubation at room temperature.
5. Wash the specimen 3 times for 2 minutes in PBS.

While rinsing, make working strength Fluorescein Antibody Solution:

Add **68ul** of blocking solution to **62 ul** of anti-Digoxigenin conjugate for a total of **130ul**. Add **370ul** of PBS to bring total volume up to **500ul**. Keep refrigerated and in the dark for *no more than 1 hour*.

6. Add Fluorescein antibody and incubate in humidifying chamber (5th floor) for 40 minutes. (Avoid exposure to light , wrap in foil and label as explained above).
7. Wash specimen 4 times in PBS, 2 minutes per wash.
8. Mount.